

DVA RAP code **AR29**

## Pressure Relieving Skin Protection



- Silicon gel provides pressure-relieving skin protection, and helps protect against skin breakage and the development of pressure sores
- Lightweight and comfortable
- Improves circulation to points on the body where pressure sores tend to form
- Uses a pressure-relieving gel found in high-end wheelchair cushions
- Cotton-lycra "coolband" stretches for ideal fit and feels great against the skin
- The coolband textile can be machine washed up to 40 deg C.
- The Gelpads can be cleaned with a damp cloth using a standard cleaning agent. They cannot be machine washed
- Supplied as pairs
- Order according to product code and size.

**Elbow & Heel**

Measure around upper arm. If for heel, measure around lower calf.

**Knee**

Measure around knee.

**Shin**

Measure around calf.

**Heel & Ankle**

Measure around ankle, just above the ankle bone.

Product Code	Description	Sizing			
		S	M	L	XL
GBHEP	Elbow & Heel	20-25	25-30	30-36	36-41
GBKP	Knee	23-28	28-33	33-41	-
GBSP	Shin	-	25-33	33-41	-
GBHAP	Heel & Ankle	20-25	25-30	30-36	-



## Heel/Elbow Protector - Foam

DVA RAP code **AR29**

- Foam padding helps maintain skin integrity and prevent ulcerations or injury to elbow or heel.
- Cushioned heel and elbow protector with foam pad can be machine washed and dried.
- Medication can be applied to pad when needed.
- Sold in pairs.

Code	Size	Description
A731S	Small	fits most small and medium women
A731M/L	Med/Large	fits larger women and average-size men.
A731XL	XLarge	fits larger men.

**Rolyan® Heel/Elbow Protector - Gel**
DVA RAP code **AR29**

- Shock-absorbing gel padding is effective following sudden injury or cumulative trauma injury to the elbow or heel.
- Helps protect the joint by diffusing and distributing pressure.
- Includes 6.4mm shock-absorbing gel pad, which won't bottom out.
- Sold in pairs.

Code	Size	Description
A730001	Small	fits most small and medium women
A730002	Med/Large	fits larger women and average-size men.
A730003	XLarge	fits larger men.

**Heel/Elbow Protector - Sheepskin**
DVA RAP code **AB11**

- The medical sheepskin heel/ elbow protector is designed to help in the healing process as well as being a preventative aid.
- Sold in pairs.
- One size.

**Application:**

- For heel and elbow pressure point control

Code	Size
CFSH05	Universal size

**DERMASAVER - Heel Elevator**
DVA RAP code **AR29**

- Elevates the heel, and covers and protects both malleoli, the bony prominences on either side of the ankle.
- Anatomically contoured and are more comfortable than other types of heel elevators.
- Soft and stretchy, with a comfort ensuring 'give'.
- DermaSaver Relevator Heel Elevators are made in 2 sizes.

Code	Size	Description
DS HES	Small	foot/heel circ below 35cm
DS HEL	Large	foot/heel circ above 35cm

