

## PULMAN Footwear



Shoe Style



Boot Style



Sandal Style



All Pulmans open as a unit, enabling ease of application.

### Features:

- Available in shoe, boot or sandal style.
- Strong enough for everyday use.
- Pulman shoes are multi-adjustable with Velcro closures, giving ease of application and a perfect, comfortable fit for feet of all shapes and sizes.
- The whole front of the shoe, plug and instep, lifts and opens as a unit. With the foot in place, the front closes and fastens securely in whatever position fits best.
- Uppers are constructed of a laminated fully lined velour fabric that stretches and breathes, maintaining optimum foot temperature.
- Pulman products have no irritating seams inside the shoes.
- The soles are light, supple, slip-resistant, adherent and quiet, with rounded edges to reduce the risk of stumbling.
- The soft foam insoles have terry textured facing, and are removable for washing or the fitting of orthopaedic insoles if required.

### Application:

- The Pulman range of therapeutic footwear is designed for those who suffer from sensitive feet or ailments which restrict mobility.
- They are particularly useful in post-operative foot care, and in the treatment of oedema, foot ulcers, rheumatoid arthritis and malformed feet or those with other foot or toe problems.
- 'Extra' fittings feature greater width and depth to accommodate swollen or bandaged feet.

### Product Care:

- Pulman footwear is sanitized for long-lasting protection against microbes and bacteria.
- Shoes and insoles are completely machine washable in warm water.
- Shoe trees are included to maintain shoe shape after washing.

### Sizing:

- Available in US shoe sizes 5 to 13.
- Order by product code and male/female US shoe size.
- Supplied in pairs.

Code	Style	Product Name	Description
NC2	Shoe	New Classic 2 Shoe	Lower fit under ankle
CX2	Shoe	Classic Extra 2 Shoe	Lower fit under ankle - with greater width and depth.
New Leiden	Shoe	New Leiden	Lower fit under ankle - with non-adjustable back.
NC1	Boot	New Classic 1 Boot	Higher fit to ankle
CX1	Boot	Classic Extra 1 Boot	Higher fit to ankle - with greater width and depth.
NL	Sandal	New Laurel	Open sandal-like design
LX	Sandal	Laurel Extra	Open sandal-like design - with greater width and depth.