



Cosmac soft foam collar

DVA RAP code AR18



Features:

- Washable polyurethane foam covered with stockinette
- Velcro fastening
- Sizing: first, measure depth (cm) from chin to clavicle for size, and then neck circumference for length
- Order by size (eg: small/regular, large/long etc)

Application:

- For treatment of mild cervical lesions
- As a reminder not to make quick, painful movements
- Ideal as night support

Product Code	Size	Depth (cm)	Fit neck circum (cm)
AF007	Small	7.5	28-35
AF011	Small/Long	7.5	35-43
AF008	Medium	8.0	28-35
AF012	Medium/Long	8.0	35-43
AF009	Large	10.0	35-43
AF014	XLarge	11.0	38-46



Philadelphia Collar

DVA RAP code AR18



Features:

- Anterior piece cradles mandible to enhance rotation control
- Made of light, strong foam, easy to trim
- Bilateral velcro fastening
- Sizing: measure depth (cm) from chin to clavicle (A, B, C, or D), and neck circumference

Application:

- For treatment of mild to moderate soft tissue damage
- Flexion and hyperextension (whiplash) injuries
- Torticollis and osteoarthritis

Depth	Neck Circumference (cm)				
	in	cm	Small (25-33)	Medium (33-40)	Large (40-46)
A	2 1/4	6	OTC6141AS	OTC6141AM	OTC6141AL
B	3 1/4	8	OTC6141BS	OTC6141BM	OTC6141BL
C	4 1/4	11	OTC6141CS	OTC6141CM	OTC6141CL
D	5 1/4	13	OTC6141DS	OTC6141DM	OTC6141DL

Order by neck size and depth (e.g. Small/A, Medium/B etc)

Exo 101 Firm Collar

DVA RAP code AR18



Features:

- Designed for symptomatic relief of neck pain or injury
- Made of polyethylene plastic, with vinyl-covered wfoam padding
- Height adjustable from 12.5 To 14.5Cm.
- Air holes for ventilation and hook-and-loop closure
- Incorporates a hinge in the bottom perimeter and bias slots under the back portion of the head.
- Automatically widens when it is lowered under the chin, and narrows when raised.
- Also available with a padded chin piece for positive positioning.

Application:

- For symptomatic relief of muscle, tendon, ligament, or soft tissue injury of the neck region.

Product Code	Size	Fit neck circum (cm)
AF101S	Small	31-37
AF101M	Medium	35-41
AF101L	Large	40-45

Sizing: measure neck circumference (cm)



Oppo 4090 Rigid Adj Collar

DVA RAP code **AR18**

Features:

- Semi-rigid support for use with soft tissue damage, osteoarthritis and whiplash injuries.
- Completely height-adjustable to provide flexion or extension positioning.
- Vinyl-covered foam padding for patients comfort and easy cleaning.

Application:

- HIVD
- Post-surgical support
- Cervical fracture
- Cervical spondylosis

To Apply:

- Adjust the height of the collar using the adjustable hook and loop straps
- Ensure the head and neck of patient are in neutral position before attaching the hook and loop closure

Product Code	Size	Fit neck circum (cm)
OPP4090S	Small	34.3-40.0
OPP4090M	Medium	40.0-46.4
OPP4090L	Large	46.4-52.7
OPP4090XL	XLarge	52.7-59.1

Measure around the circumference of the neck.



Oppo 4190 Rigid Adj Collar w/ chin

DVA RAP code **AR18**

Features:

- Completely height-adjustable to provide desired flexion or extension positioning.
- Vinyl-covered foam padding for patient's comfort and easy cleaning.
- Ventilated to allow air circulation.
- Hook and loop closure for quick, easy application and removal.

Application:

- HIVD
- Post-surgical support
- Cervical fracture
- Cervical spondylosis

To Apply:

- Apply the collar to the neck with the chin on the chin rest.
- Adjust the height of the cervical collar using the hook and loop straps to decide the proper position of the head.
- Before attaching the hook and loop closure strap, make sure the head and neck are in neutral position.
- Fasten the hook and loop strap securely.

Product Code	Size	Fit neck circum (cm)
OPP4190S	Small	34.3-40.0
OPP4190M	Medium	40.0-46.4
OPP4190L	Large	46.4-52.7

Measure around the circumference of the neck.




Features:

- Designed to optimize support and comfort, the key components for better patient outcomes.
- Provides substantial motion restriction without producing painful pressure points that can lead to skin breakdown or poor patient compliance.
- All contact surfaces of the collar are cushioned with cotton-lined, breathable foam padding.
- For detailed information regarding sizing, fitting, replacement pads & accessories, visit www.aspenmp.com

Application:

- Sprain/strain
- Radiculopathy
- Disc syndrome
- Musculoskeletal pain
- Trauma / whiplash
- Fracture management
- Cervicogenic headache
- Cervical disc syndrome
- Decompressive procedures
- Post surgical stabilization

DVA RAP code **AR18**


Collar	Collar Set	Replacement pads code	Size	Height (cm)	Neck Circum (cm)
(1 set of pads)	(2 sets of pads)				
AS983108	AS983128	AS983159	Adult short	2	33-53
AS983110	AS983130	AS983159	Adult regular	2-3	33-53
AS983112	AS983132	AS983159	Adult tall	3-5	33-53
AS983114	AS983134	AS983159	Adult X/tall	5+	33-53

Height measurement: chin to shoulder as per illustration



Miami-J Collar

DVA RAP code AR18



Features:

- A phenotype-driven sizing system — with unique fronts and backs that are bioengineered to work together for superior immobilization — makes size selection simple, while drastically reducing the margin of error in fitting a patient
- Adjustable mandibular and occipital support enables customization to disburse pressure and avoid skin breakdown
- Sorbatex™ pads that are specifically engineered to be anti-decubitus — protecting skin integrity during extended wear — and are removable for easy cleaning
- Sternal pad disburses pressure and is removable for comfort and improved compliance when supine and/or swallowing

- Patient assist strap allows for flexibility and comfort without sacrificing immobilization
- Latex-free and compatible with all imaging modalities
- For detailed information regarding sizing, fitting, replacement pads & accessories, visit www.Ossur.Com

Application:

- C-spine precaution for trauma patients
- Immobilization for pre and post c-spine surgery
- Degenerative disorders
- Spinal stenosis
- Spondylolisthesis

Collar (1 set of pads)	Collar Set (2 sets of pads)	Replacement pads code	Size	Height (cm)	Neck circ (cm)	Sizing Guidelines
Standard Sizes: > 70% of Population						
MJ-300	MJR-300	MJP-100	Small	5.0	25-50	<ul style="list-style-type: none"> • standard for mature men • short-necked women
MJ-400	MJR-400	MJP-100	Medium	6.3	25-50	<ul style="list-style-type: none"> • standard for mature women • thin, mature men
MJ-500	MJR-500	MJP-100	Large	7.5	25-50	<ul style="list-style-type: none"> • long tall neck • 'swan' neck • young women • adolescents
Specialty Sizes: < 30% of Population						
MJ-200L	MJR-200L	MJP-200L	Stout	4.4	35-60	<ul style="list-style-type: none"> • very large neck circumference • 'no neck' • obese • massive shoulders
MJ-200S	MJR-200S	MJP-100	Super Short	2.5	25-50	<ul style="list-style-type: none"> • kyphotic (chin-on-chest) neck • osteoporotic and/or ankylosing spondylitic • geriatric
MJ-250	MJR-250	MJP-250	XSmall	3.1	20-35	<ul style="list-style-type: none"> • very short, thin neck principally female • more prevalent in Asian population

Measurements: height of neck - chin to shoulder (as per Aspen collar)



Body Assist N71 Thermal Shoulder Brace

DVA RAP code AR03

Neoprene (black)



Features:


- Provides soft tissue support & compression.
- Increases elasticity of muscles, ligaments & tendons.
- Assists in prevention of injury or re-injury.
- Elastic support strap fits comfortably under the arm pit & assists to generate better than average compression & support to the shoulder.
- Multi-directional compression assists to reduce soft tissue swelling.
- Can be worn while playing sports.
- Promotes greater blood circulation, provides a more speedy recovery.
- Contoured panels fit most body shapes.
- Completely washable & hygienic.

Stabilizer strap:

- The stabilizer strap is a rigid woven strap that adheres with hook/loop over the top of the body assist n71b shoulder brace.
- It can be applied in various configurations to minimize shoulder movement and enhance stability.

Application:

- Adjustable support for symptomatic relief of sprains, strains, bursitis & tendonitis.
- May be used alone, or with optional stabilising strap for added support (product code n71ac)



STABILIZER STRAP

Product Code	Size	To Fit
NHN71AC1	Regular	To fit N71B brace XS to M
NHN71AC3	Large	To fit N71B brace M to XL

SHOULDER BRACE

Left	Right	Size	cm
NHN71L0	NHN71R0	x small	80 - 90
NHN71L1	NHN71R1	small	90 - 100
NHN71L2	NHN71R2	medium	100 - 110
NHN71L3	NHN71R3	large	110 - 120
NHN71L4	NHN71R4	x large	120 - 130

Measure circumference around chest, around outside the arm on affected side.




Rehband 7126 Shoulder Support

Neoprene (black)

DVA RAP code

AR03



Features:

- Soft and flexible shoulder support with excellent fit.
- Provides warmth and pain relief to the shoulder joint.

Application:

- Used for over-loaded or degenerative shoulders, inflammations

Left	Right	Size	cm
REH7126XS-L	REH7126XS-R	x small	24-28
REH7126S-L	REH7126S-R	small	28-32
REH7126M-L	REH7126M-R	medium	32-35
REH7126L-L	REH7126L-R	large	35-38
REH7126XL-L	REH7126XL-R	x large	38-42

Sizing: measure biceps circ (cm)




Oppo 1072 Shoulder Support

Neoprene (beige)

DVA RAP code

AR03



Features:

- Provides compression and support to shoulder and aids healing.
- Therapeutic warmth for the shoulder and upper arm.
- Offers support and protection.
- Minimizes potential for future injuries
- Fits left or right.

Application:

- Rotator cuff tendinitis
- Bursitis
- Shoulder instability

Product Code	Size	cm.
OPP1072S	S	83.2 - 93.3
OPP1072M	M	93.3 - 103.5
OPP1072L	L	103.5 - 113.7
OPP1072XL	XL	113.7 - 126.4

Measure around the chest (cm)




Oppo 2172 Reinforced Shoulder Support

Neoprene (beige)

DVA RAP code

AR03



Features:

- Two reinforced elastic straps on upper arm provide extra support for the shoulder.
- Allows the patient gradually return to full range of motion
- Fits left or right.

Application:

- Rotator cuff tendinitis
- Bursitis
- Shoulder instability

Product Code	Size	cm.
OPP2172S	S	81.6 - 91.8
OPP2172M	M	91.8 - 101.9
OPP2172L	L	101.9 - 114.6
OPP2172XL	XL	114.6 - 127.3


Measure around the chest (cm)

Barrere Clavicle Brace

(white)

DVA RAP code

AR03



Features:

- Padded 'figure of eight' shoulder straps
- One size fits all

Application:

- Provides control of forward or downward displacement of the shoulders

Product Code	Size
B422	Universal Size


One size fits all



Oppo 2075 Clavicle Brace

(beige)

DVA RAP code AR03



Features:

- Excellent aid in correcting poor posture and shoulder slump.
- Adjustable hook and loop closure facilitates proper fit.

Application:

- Poor posture
- Weak scapula muscles

To Apply:

- Bring the hook and loop strap through the metal ring to form a loop on each side. Slip the arm through the loop. Keeping the straps in the back.

- With padded loops over the shoulders, pull both straps to the front.
- Fasten straps snugly.

Product Code	Size	cm.
OPP2075S	S	53.7-63.8
OPP2075M	M	63.8-76.5
OPP2075L	L	76.5-89.2
OPP2075XL	XL	89.2-104.5
OPP2075XXL	XXL	104.5-117.2
OPP2075XXXL	XXXL	117.2-129.9

Measure around scapula under arm.



Oppo 4072 humerus brace

(beige)

DVA RAP code AR03



Features:

- Padded webbing strap fits under the opposite axilla
- Hook and loop closure for quick, easy, adjustable application and removal.

Application:

- Stabilizes and supports the injured shoulder

Product Code	Size	cm.
OPP4072S	S	24-29
OPP4072M	M	29-34
OPP4072L	L	34-39

Measure mid upper arm circumference (cm)



Rolyan hemi arm sling II

(white)

DVA RAP code AR03



Features:

- Approximately 3.8cm longer than the original model
- Effective unilateral anti-subluxation support.
- Positions the head of the humerus in the glenoid fossa and the scapula in external rotation.
- Enhances normal gait pattern.
- Elastic straps are dynamic, allowing a patient to actively move returning muscles and perform self-ranging exercises.
- The supports enable hemiplegics to use their affected arm as an assist during functional tasks.
- Easy to don/doff and can be worn under clothing.
- Optional universal cuff adapter A653-11 available - positions the hand above heart to treat or prevent edema. Attaches to a d-ring on the unaffected side.

Application:

- Effective following stroke, soft tissue injury or ligament strain.
- Provides support for heavier patients as well as patients with increasing active motion"

Left	Right	Size	cm	chest
A6531	A6532	small	36-46	-
A6533	A6534	medium	46-61	76-89
A6535	A6536	medium long	46-61	89-102
A6537	A6538	large	61-71	-
A6539	A65310	x large	71-86	-

To fit, measure the distance from acromion to acromion posteriorly across the midpoint of the scapula. To determine if patient is medium or medium long, measure the chest circumference at the axilla.



Rolyan MFC II shoulder orthosis

DVA RAP code AR03



Features:

- lycra/spandex fabric offers excellent conformability to body contours for a snug fit
- easy to put on and take off
- can be worn under clothing
- hand wash, air dry

Application:

- provides support for weak shoulder muscles and strained tendons or ligaments
- limits shoulder motion and maintains proper humerus positioning to allow healing

Left	Right	Size	cm
A54551	A54552	S	23-28
A54561	A54562	M	28-33
A54571	A54572	L	33-38

Measure circumference of mid upper-arm.




Rehband 5055 Acro Comfort Shoulder Support

Neoprene (black)

DVA RAP code

AR03



Features:

- Stabilises the shoulder joint.
- Stretch flex system is anatomically designed for optimal function.
- Excellent comfort through the temperature regulating features of comfortemp®.
- Fits left or right

Application:

- Instability and hyper mobility.

Product Code	Size	cm
REH5055XS	x small	24-28
REH5055S	small	28-32
REH5055M	medium	32-35
REH5055L	large	35-38
REH5055XL	x large	38-42

Measure circumference around the shoulder joint going under the arm. (cm)



Rehband 1068 EQ Posture Support

Elastic (grey)

DVA RAP code

AR03



Features:

- A flexible aid for training a better posture
- Bendable stays for additional support
- Elastic and well-fitting
- Wear with singlet or t-shirt underneath

Application:

- Used for training when the posture is poor or the thoracic muscles are weak

Product Code	Size	cm
REH1068XS	x small	56 - 64
REH1068S	small	64 - 72
REH1068M	medium	72 - 80
REH1068L	large	80 - 90
REH1068XL	x large	90 - 100

Sizing: waist circ (cm)




Oppo 2175 Posture Aid

Elastic (beige)

DVA RAP code

AR03



Features:

- Triangle configuration correctly distributes pressure to tip of shoulder.
- Suitable for adults or children.
- Adjustable hook and loop closure facilitates proper fit.

Application:

- Poor posture
- Weak scapula muscle

Product Code	Size	cm.
OPP2175S	S	64.1-76.5
OPP2175M	M	76.5-89.2
OPP2175L	L	89.2-104.5
OPP2175XL	XL	104.5-117.2
OPP2175XXL	XXL	117.2-129.9

Measure around scapula under arm.