



TRIBUTE™
MANUFACTURED BY SOLARIS

Custom Nightwear Lymphoedema Garments

Why is Tribute essential?

- Tribute custom-designed NightWear garments ensure a successful lymphoedema home management program.
- Tribute garments greatly simplify the patient's home programme, decrease treatment time, and increase patient compliance.
- Worn nightly, Tribute garments eliminate the need for complex bandaging, and offer an effective solution to help regulate fluid fluctuations, decrease tissue turgor, assist a weakened lymphatic system, and maintain gains made in therapy.

How does Tribute work?

- Foam-filled channels in the garments create localised tissue stretch and pressure differential, opening the initial lymphatic gaps to encourage re-absorption of interstitial fluids and particles.
- Every Tribute garment is designed to maximise therapeutic results: Vertical style channels provide directional flow, distal to proximal; Chevron style provides bi-directional flow – medial to lateral and distal to proximal, directing interstitial fluid towards collateral pathways.

Comfort:

- The proprietary fabric used in every Tribute garment, ensures tissue health by providing moisture and bacterial management, and thermo release properties to keep the wearer cool and comfortable.
- Unlike bandaging and other nightwear garments, Tribute allows the wearer the comfort of a full range of motion, allowing sleep in any position.

Contraindications:

- Arterial insufficiency or degeneration.
- deep vein thrombosis.
- untreated congestive heart failure, localized systemic infection or cancer.
- absent or severely impaired sensation.


Warranties:

- Sure-Fit Warranty: within two weeks of receipt, if the garment does not fit as desired, Solaris will alter or remake the garment free of charge.
- Alteration Policy: Solaris offers one free alteration to decrease the size of the garment following girth reduction. Unfortunately, it is not possible to enlarge the garment.
- One Year Warranty: Solaris extends a one-year warranty against any defect in workmanship and material.

Further Information:

- Due to the extensive range of Tribute garments, a separate detailed Tribute catalogue, including measuring charts, is available from Cosmac Healthcare on request.
- Below is a cross-section of custom Tribute nightwear lymphoedema garments and garment combinations.

TRIBUTE Upper Extremity



Glove



Extended Glove



F/Tips to Axilla



F/Tips to Clavicle



F/Tips to Clavicle with Vest



Wrist

MCPs

F/tips

Sleeves may be from the wrist, MCPs or F/tips, to the axilla or clavicle. Wrist to axilla/glove combinations may also be used

TRIBUTE Lower Extremity



Toes to Knee




Toes to Mid-Thigh




Toes to Groin



Chap and Half-Leg



Pant



Pant with Foot

TRIBUTE Trunk & Torso



Shorts



Lower Torso



Abdominal Binder



Bustier



Vest



Body Suit



Thoracic Sling Wrap



Short Sleeve Shirt



Upper Torso/Neck


TRIBUTE Head & Neck




Cervical Collar



Mandibular Unit



Orbital Unit



Mandibular/Orbital Unit

BREAKING NEWS

As this catalogue goes to print, we are pleased to announce the following exciting SOLARIS product additions, which will come on stream in early 2013.

- **TributeDay** – this new custom-made flat-knit garment for day-time wear, will complement the existing Tribute nightwear garment (which will then be aptly renamed TributeNight!). The addition of this product will then provide a made-to-measure 24 hour solution for lymphoedema management.

- **ExoStrong** – off-the-shelf flat-knit hand, arm & leg garments. Available in 20-30 and 30-40mmHg, these will offer a more economical solution than made-to-measure garments for those with standard measurements.