



### Oppo 1038 knee support – Wrap around → Open patella → Stabilising flexible stays

Neoprene (black)

DVA RAP code AR02



#### Features:

- two stays on both sides provide additional stability
- adjustable hook and loop closure for custom fit
- provides warmth and equal compression to weak or overstress knee

#### Application:

- slight instability of the knee
- articular effusions and swellings
- mild strains of the knee
- preventative use in sports

#### To Apply:

- loosen straps and wrap the support around the knee joint
- Position the knee cap in the patella opening
- fasten the hook and loop closures above and below the patella snugly

Product Code	Size	cm
OPP1038S/M	S/M	33.0-38.1
OPP1038L/XL	L/XL	38.1-50.8

Measure around the mid-section of the knee, while in a standing position



### Oppo 1230 adjustable knee stabilizer – Wrap around → Open patella → Stabilising flexible stays

Coolprene (beige)

DVA RAP code AR02



#### Features:

- two medial and lateral spiral stays provide additional support and stability for the knee
- the additional circle pad around the patella provides padded protection and stabilization to prevent patella displacement
- open patella design to relieve pressure
- heat retention and compression reduces potential for strain or injury

#### Application:

- knee instability
- mild strain/ sprain
- collateral ligament injury

#### To Apply:

- position the opening over kneecap and be sure the two lateral stays are on each side of knee
- wrap and secure the shorter strap behind knee joint, between upper and lower straps
- wrap and secure the upper and lower straps to front of support

OPP1230	One size	31.8-40.6
---------	----------	-----------

Measure around the knee joint.



### Oppo 1130 knee stabiliser adjustable – Wrap around → Open patella → Stabilising flexible stays

Neoprene (beige)

DVA RAP code AR02



#### Features:

- two medial and lateral spiral stays provide additional support and stability for the knee
- heat retention and compression reduces potential for strain or injury
- adjustable for custom fit

#### Application:

- knee instability
- mild strain/ sprain
- collateral ligament injury

#### To Apply:

- position the opening over kneecap and be sure the two lateral stays are on each side of knee
- wrap and secure the shorter strap behind knee joint, between upper and lower straps
- wrap and secure the upper and lower straps to front of support

OPP1130	One size	31.8-40.6
---------	----------	-----------

Measure around the knee joint.



## Rehband 8362N genu carezza knee support – Wrap around → Open patella → Stabilising flexible stays

Elastic (grey)

DVA RAP code AR02



### Features:

- support and stabilization by coil springs
- supports sensory motor control
- wraparound design for easy application
- lateral and medial coil springs for stability and to prevent the brace from rolling
- anatomic design with soft trimmings to provide excellent comfort and fit
- innovative material that combines compression and breathability
- two elastic straps
- can be used for left or right

### Application:

- knee pain with no or minimal ligament laxities and no or light perception of instability
- soft tissue irritations (e. g. tendomyopathies, insertion tendinitis, swellings, light effusions)

- overuse (e. g. sports)
- posttraumatic and postsurgical irritations
- beginning osteoarthritis
- patellofemoral pain syndrome (e. g. chondromalacia, chondropathy of the patella, beginning osteoarthritis of the patellofemoral joint)

Product Code	Size	15 cm below mid patella	15cm above mid-patella
REH8362NXXS	XXS	29-32	36-40
REH8362NXS	XS	32-35	40-44
REH8362NS	S	35-38	44-48
REH8362NM	M	38-41	48-52
REH8362NL	L	41-44	52-56
REH8362NXL	XL	44-48	56-61
REH8362NXXL	XXL	48-51	61-67

Measure around leg, 15cm above, and 15cm below middle of knee.



## Body Assist 43C elastic cartilage knee brace – Wrap around → Open → Stabilising flexible stays

Elastic (beige)

DVA RAP code AR02



### Features:

- adjustable front closure, lateral/medial condyle pads & patella cutout.
- firm elastic orthopaedic soft tissue support with additional compression over the medial and lateral meniscus (cartilage).
- three medial and three lateral spring steel stays provide greater support.

### Application:

- supports the medial and lateral ligaments and relieves pressure from the patella.

Product Code	Size	cm
NH43C1	S	To 32
NH43C2	M	32-37
NH43C3	L	37-42
NH43C4	XL	42-47
NH43C5	XXL	47 +

Measure around the knee joint.



## Rehband 8388 genu therma stable knee support – Pull on → Open Patella → Stabilising flexible stays

Neoprene (grey)

DVA RAP code AR02



### Features:

- stabilizing knee brace with a pad providing extra support for the kneecap
- knee support made from warming and flexible neoprene
- medial and lateral springs provide stabilization
- can be used for left or right

### Application:

- knee pain without ligament laxities or perception of instability
- soft tissue irritations (e. g. tendomyopathies, insertion tendinitis, swellings, light effusions)
- overuse (e. g. sports)
- long-term soft tissue irritations after injuries and / or surgery

- beginning osteoarthritis
- patellofemoral pain syndrome (e.g. chondromalacia, beginning osteoarthritis of the patellofemoral joint)

Left	Size	15 cm below mid patella	15 cm above mid patella
REH8388XXS	XXS	29-32	36-40
REH8388XS	XS	32-35	40-44
REH8388S	S	35-38	44-48
REH8388M	M	38-41	48-52
REH8388L	L	41-44	52-56
REH8388XL	XL	44-48	56-61
REH8388XXL	XXL	48-51	61-67

Measure around leg, 15cm above, and 15cm below middle of knee.



### Oppo 1030 knee stabilizer – Pull on → Open patella → Stabilising flexible stays

Neoprene (beige)



#### Features:

- neoprene circle pad around the patella provides padded protection and stabilization to prevent patella displacement.
- two medial and lateral stays provide additional support and stability for the knee.
- two hook and loop straps insure comfortable fit and stabilization during activities.
- helps relieve pain from tendinitis and arthritis.
- retains body heat to warm knee joint and increases blood circulation.

#### Application:

- arthritis
- knee strain/ sprain
- knee instability
- patellar instability

DVA RAP code AR02

#### To Apply:

- hold the top of support using both hands, one on each side.
- insert foot into support and pull up to the knee.
- position the opening over the kneecap.
- secure the upper small hook and loop strap.
- be sure the two medial and lateral stays are on each side of knee.
- thread through buckles and secure the upper and lower hook and loop straps to desired compression.

Product Code	Size	cm
OPP1030S	S	33.7-37.5
OPP1030M	M	37.5-40.6
OPP1030L	L	40.6-45.1
OPP1030XL	XL	45.1-50.2

Measure around the mid-section of the knee, while in a standing position.



### McDavid A415B cartilage knee support – Pull on → Open patella → Stabilising flexible stays

Neoprene (black)



#### Features:

- buttress protects and relieves pressure from the patella.
- two medial and 2 lateral spring steel stays provide greater support.
- firm compression and soft tissue support with additional compression over the medial and lateral ligaments.
- cushions against minor scrapes and bumps.

#### Application:

- supports the medial and lateral ligaments and relieves pressure from the patella.

DVA RAP code AR02

Product Code	Size	cm
NHA415B1	S	30-33
NHA415B2	M	33-38
NHA415B3	L	38-43
NHA415B4	XL	43-51
NHA415B5	XXL	51-56

Measure around the knee joint.



### McDavid A416 thermal collateral knee support – Pull on → Open patella → Stabilising flexible stays

Neoprene (black)



#### Features:

- sewn-in buttress protects and relieves pressure from the patella.
- includes 2 medical spring stays on both sides of the knee for extra support.

#### Application:

- supports the medial and lateral ligaments and relieves pressure from the patella.

DVA RAP code AR02

Product Code	Size	cm
NHA4161	S	30-33
NHA4162	M	33-38
NHA4163	L	38-43
NHA4164	XL	43-51
NHA4165	XXL	51-56

Measure around the knee joint.



**McDavid A421 patella knee support** – Pull on → Open patella → Stabilising flexible stays  
Neoprene (black, beige)

**DVA RAP code** AR02



**Features:**

- delivers firm compression & soft tissue support.
- made of medical grade lightweight neoprene that retains body warmth.
- contoured design for a better fit with two medial and two lateral medical steel stays.
- top and bottom anchor straps to minimize slippage.
- available in Beige or Black

Black	Beige	Size	cm
NHA421B1	NHA421S1	S	30-33
NHA421B2	NHA421S2	M	33-38
NHA421B3	NHA421S3	L	38-43
NHA421B4	NHA421S4	XL	43-51
NHA421B5	NHA421S5	XXL	51-56

Measure around the knee joint.



**McDavid A425 ligament knee support** – Pull on → Open patella → Stabilising flexible stays  
Neoprene (black, beige)

**DVA RAP code** AR02



**Features:**

- delivering firm compression and soft tissue support.
- made of medical grade lightweight neoprene that retains body warmth.
- 4 elastic straps that deliver additional support to ligaments.
- horse-shoe buttress to stabilise the knee cap.
- contoured design for a better fit with two medial and two lateral medical steel stays.
- top and bottom anchor straps minimize slippage

Black	Beige	Size	cm
NHA425B1	NHA425S1	S	30-33
NHA425B2	NHA425S2	M	33-38
NHA425B3	NHA425S3	L	38-43
NHA425B4	NHA425S4	XL	43-51
NHA425B5	NHA425S5	XXL	51-56

Measure around the knee joint.



**Oppo 2233 spiral knee brace** – Pull on → Open patella → Stabilising flexible stays  
Elastic (beige)

**DVA RAP code** AR02



**Features:**

- promotes healing, keeps joints and muscles flexible
- two medial and lateral spiral stays on both sides provide additional support and stability for the knee
- patella opening eliminates pressure on kneecap
- helps relieve pain resulting from injured or weak knee conditions

**Application:**

- weak knee
- knee strain/ sprain

**To Apply:**

- hold the top of brace using both hands, one on each side
- insert foot into brace and pull up to the knee
- the opening of brace should be over kneecap
- each two stays should be on each side of knee

Product Code	Size	cm
OPP2233S	S	30.5–34.3
OPP2233M	M	34.3–38.1
OPP2233L	L	38.1–41.9
OPP2233XL	XL	41.9 – 45.7

Measure around the mid-section of the knee, while in a standing position



**Body Assist 42S knee support** – Pull on → Open patella → Stabilising flexible stays  
Elastic (beige)

**DVA RAP code** AR02



**Features:**

- an open patella with medial & lateral medical side stays for extra overall support
- open Patella minimises pressure over the knee cap
- quality lightweight elastic orthopaedic soft tissue support with additional side supports

**Application:**

- weak knee
- knee instability
- knee strain/ sprain

Product Code	Size	cm
NH42S1	S	30.5-35.5
NH42S2	M	33-38
NH42S3	L	35.5-40.5
NH42S4	XL	40.5-46
NH42S5	XXL	46 +

Measure around the knee joint.

**BodyAssist**
**Body Assist 41S cross-cut elastic knee brace** – Pull on → Closed patella → Stabilising flexible stays  
 Elastic (beige)
**Features:**

- anatomically shaped by a cross-cut design that minimizes irritation & pressure behind the knee.
- provides orthopaedic support for minor knee joint strains & sprains.
- removable medical support spring stays on sides of the knee joint to provide medio-lateral support.
- brace retains some body warmth & can increase blood circulation.
- lightweight, ventilated elastic for all day comfort.

**Application:**

- weak knee
- knee instability
- knee strain/ sprain

DVA RAP code AR02

Product Code	Size	cm
NH41S1	S	To 32
NH41S2	M	32-37
NH41S3	L	37-42
NH41S4	XL	42-47
NH41S5	XXL	47 +

Measure around the knee joint.