



## Dale Abdominal Binder

DVA RAP code AR08



### Features:

- single ply stretch loop elastic (latex-free)
- velcro closure anywhere on belt for perfect fit
- provides full, even pressure and aids abdominal muscle
- colour white
- Sizing: available in either 23cm (pictured) or 30cm widths

### Application:

- suitable for post-surgical / incisional hernia
- general abdominal support

### Product care:

- join the hook and loop fasteners before washing
- hand wash in warm water using mild detergent
- dry out of direct sunlight

23 cm wide	30cm wide	Size	cm
DO 410	DO 810	Regular	to 114
DO 411	DO 811	Large	114 - 157
DO 418	DO 818	X Large	152 - 191

Sizing: measure circumference of waist.



## Scudotex 644 Elastic hernia underpants w/pads

DVA RAP code AR14



### Features:

- designed with special highly compressive fabric
- anatomical shape ensures support for the abdominal area to prevent muscle prolapse and provide a balanced inguinal compression
- foam pads fit into pockets to support inguinal hernia

### Application:

- comfortable and effective in cases of single or double inguinal hernia
- for pre or post-operative use

Product Code	Size	waist	hip
SCTX644 3	M	80-86	90-98
SCTX644 4	L	86-92	96-104
SCTX644 5	XL	92-99	102-111
SCTX644 6	XXL	99-106	109-118
SCTX644 7	XXXL	106-114	116-126

Measure hip and waist circumference (cm)



## Scudotex 609 Hernia support belt w/pads

DVA RAP code AR14



### Features:

- inguinal truss with anatomic protective front in hygienic elastic and cotton fabric
- adjustable leg straps go under the thighs and attach at the back to allow correct inguinal compression
- designed with two removable rubber pads covered with hygienic cotton
- back adjustable opening

### Application:

- comfortable and effective in cases of single or double inguinal hernia
- for pre or post-operative use
- in cases of post operative recovery it is recommended to remove both pads

Product Code	size	cm
SCTX609 2	S	75-80
SCTX609 3	M	80-86
SCTX609 4	L	86-92
SCTX609 5	XL	92-99
SCTX609 6	XXL	99-106
SCTX609 7	XXXL	106-114
SCTX609 8	XXXXL	114-122

Measure waist circumference (cm)



## Body Assist 590 hernia belt

DVA RAP code AR14



### Features:

- removable foam pads
- under leg straps
- easy velcro fastening and adjustment
- made of lightweight material construction & assists to inhibit skin irritation
- the belt can be worn whilst showering or swimming

### Application:

- comfortable and effective in cases of single or double inguinal hernia
- for pre or post-operative use.
- in cases of post operative recovery it is recommended to remove both pads

### To Apply:

- in case your hernia is one-sided, you will need to remove the pad not used. However, many doctors recommend keeping both pads in place so that even pressure is applied to both sides
- position yourself lying down on your back
- gently massage the swelling until it retracts back inside the abdominal wall

- now place the pressure pads over the affected area
- do not remove the pads from their pockets
- secure the waist strap so that the fit is firm
- secure the leg straps and stand up
- now do final adjustments for maximum comfort
- make sure the pad(s) stay over the affected area
- it is advisable to wear the belt over an undergarment

### Product care:

- washing instructions: remove pads and join the hook and loop fasteners
- hand wash in warm water using mild detergent
- do not bleach, tumble dry or dry-clean
- do not dry in direct sun-light

Product code	Size	cm
NH590 0	x small	65 to 75
NH590 1	small	75 - 85
NH590 2	medium	85 - 95
NH590 3	large	95 - 105
NH590 4	x large	105 - 115

Sizing: measure waist circumference (cm)



## Scudotex 682 Inguinal hernia belt – heavy duty

DVA RAP code AR14



### Features:

- quality inguinal truss in strong elastic
- back adjustable
- soft cotton edges for greater comfort
- leather-covered pads
- adjustable leg straps
- double pads

### Application:


- heavy duty hernia belt for all inguinal hernias

Product Code	Size	waist
SCTX682 2	S	75-80
SCTX682 3	M	80-86
SCTX682 4	L	86-92
SCTX682 5	XL	92-99
SCTX682 6	XXL	99-106

Measure waist circumference (cm)

## Barrere 24 post-op inguinal belt

DVA RAP code AR14



### Features:

- heavy duty post-operative inguinal hernia support belt
- same as model 25-S, but with no pads fitted
- therefore, no left/right fitting
- two understraps ensure correct positioning

### Application:

- post-operative inguinal hernia support belt

### Measuring:

- take body circumference 5cm below the hips, with tape measure sloping down at the front towards the abdomen

Product Code	Size	cm
B24 XS	x small	76 - 81
B24 S	small	81 - 86
B24 M	medium	86 - 91
B24 L	large	91 - 96
B24 XL	x large	96 - 101
B24 XXL	xx large	101 - 106

Sizing: measure circumference (cm) 5cm below the hips, with tape measure sloping down at the front towards the abdomen.

## Barrere 25-S inguinal hernia belt

DVA RAP code AR14



### Features:

- heavy duty hernia belt for all inguinal hernias
- available in left, right or double
- two understraps ensure correct positioning

### Application:

- heavy duty inguinal hernia support.

### Measuring:

- take body circumference 5cm below the hips, with tape measure sloping down at the front towards the abdomen

Product Code	Size	cm
B25-S XS	x small	76 - 81
B25-S S	small	81 - 86
B25-S M	medium	86 - 91
B25-S L	large	91 - 96
B25-S XL	x large	96 - 101
B24-S XXL	xx large	101 - 106

Sizing: measure circumference (cm) 5cm below the hips, with tape measure sloping down at the front towards the abdomen.

## Barrere 26-S scrotal hernia belt

DVA RAP code AR09



### Features:

- available in left, right or double
- two understraps ensure correct positioning
- with suspensory and double-adjustable crossover straps to ensure complete control of any size reducible hernia

### Application:

- heavy duty hernia belt for all scrotal hernias

### Measuring:

- take body circumference 5cm below the hips, with tape measure sloping down at the front towards the abdomen

Product Code	Size	cm
B26-S XS	x small	76 - 81
B26-S S	small	81 - 86
B26-S M	medium	86 - 91
B26-S L	large	91 - 96
B26-S XL	x large	96 - 101
B26-S XXL	xx large	101 - 106

Sizing: measure circumference (cm) 5cm below the hips, with tape measure sloping down at the front towards the abdomen.



## McDavid A310 athletic support

DVA RAP code AR09



### Features:


- a multi-purpose supporter featuring a power-knit pouch with Intera for rapid drying and cooling
- provides greater resistance to bacteria and mildew, and the odours and skin irritation they produce
- more comfortable than cotton or synthetics
- Incorporates an 8cm knit waistband and 1cm leg straps
- high rubber content assures maximum support and durability

Product Code	Size	cm
NHA310 1	small	66-82
NHA310 2	medium	82-97
NHA310 3	large	97-112
NHA310 4	x large	112-128

Measure around waist (cm)

**McDavid** McDavid A321 athletic support with protective cup

DVA RAP code AR09



**Features:**

- a durable and comfortable supporter featuring a plastic protector cup, enclosed in a power-knit pouch, with Intera for rapid drying and cooling
- provides greater resistance to bacteria and mildew, and the odours and skin irritation they produce
- more comfortable than cotton or synthetics
- Incorporates an 8cm knit waistband and 1cm leg straps
- high rubber content assures maximum support and durability

Product Code	Size	cm
NHA321 1	small	66-82
NHA321 2	medium	82-97
NHA321 3	large	97-112
NHA321 4	x large	112-128

Measure around waist (cm)

**Supreme Suspensory lightweight**

DVA RAP code AR09



**Features:**

- light-weight scrotal support
- designed to hold and lift the scrotum
- relieves the strain and fatigue of swollen or sagging testicles
- opening for penis

**Application:**

- for use post-surgical, for injury, or to support the cremaster muscle

Product Code	Size	cm
ASO 13	small	76 - 89
ASO 14	medium	89 - 104
ASO 15	large	104 - 117
ASO 15XL	x large	117 - 130

Measure around waist (cm)

**TRUFORM®** Truform 243 suspensory - heavy duty

DVA RAP code AR09



**Features:**

- for larger scrotal swelling/enlargement
- heavy duty construction
- 4.5cm elastic web waistband with claw-buckle fastening.
- strong mesh pouch
- 2cm elastic leg-straps
- opening for penis
- sized according to pouch
- medium is the size most often selected

**Application:**

- relieves the strain and fatigue of swollen or sagging testicles

Product Code	Size
TR243M	medium
TR243L	large
TR243XL	x large

**TSL Scrotal Support – extra large**

DVA RAP code AR09



**Features:**

- for extreme scrotal swelling/enlargement provides support and patient comfort by alleviating heavy sensation
- comfortable to wear
- soft Powernet fabric

**Application:**

- scrotal oedema
- scrotal lymphoedema
- scrotal hernia
- scrotal rupture

**Sizing:**

- contact Cosmac to determine if patient can be fitted with standard sizes, or requires a custom made garment

Product Code	Size	Combined Scrotum Size
TSLSS 1	small	orange
TSLSS 2	medium	rock melon
TSLSS 3	large	watermelon