



## Rolyan palm protector with finger separator

**DVA RAP code** AR03


### Features:

- made from closed-cell foam, lined with pile
- soft finger separator strap weaves between fingers and can be cut to adjust fit
- hand washable
- one size fits all

### Application:

- prevents fingers from digging into palm while abducting fingers to decrease tone

Left	Right
A81201	A81202

Universal size, left and right.



## Rolyan palm shield - foam

**DVA RAP code** AR03


### Features:

- this foam support wraps around the user's wrist and hand and holds the thumb and fingers in place
- it is held in place with a Velcro closure at the wrist
- the support has four rounded foam sections that separate and position the fingers
- it also has a padded section that rests in the user's palm, and supports and positions the user's thumb

- this hand support is designed for users with moderate to severe hand tightness
- One size fit most adults – left and right fittings

Left	Right
A812401	A812402

Universal size, left and right.



## Rolyan palm shield - sof-gel

**DVA RAP code** AR03


### Features:

- a palm shield with gel which helps prevent skin breakdown in the palm and between the fingers
- gel in finger separators and palm won't bottom out giving the product a longer life
- ideal for moderate to severe hand tightness
- helps prevent finger contractures
- large thumb piece provides excellent thumb positioning
- One size fits most adults

Left	Right
A812405	A812406

Universal size, left and right.



## Rolyan palm protector

**DVA RAP code** AR03


### Features:

- holds fingers out to prevent digging into the palm of the hand
- cuff fabricated from closed-cell foam lined with a pile fabric
- pre-rolled and attached to simplify application
- hand wash, air dry
- one size for each of left and right hand

Left	Right
A812101	A812201

Universal size, left and right.

## Rolyan® Rolyan Carve-It hand positioning splint

DVA RAP code AR03



### Features:

- for use with strokes, Alzheimer's disease, dementia and other neurological impairments
- helps anatomically position the hand with mild to moderate finger spasticity
- design allows use with left and right hand
- easy to modify by cutting foam with an electric knife for a custom fit
- foam shaped to fit the hand and comfortable for the patient
- unlike traditional hand cones, this splint maintains the

- arches of the hand in an appropriate position
- flat bottom provides easy positioning in bed or on lap tray, arm tray, etc
- removable cover can be hand washed and air dried

Splint	Replacement Cover	Size	Cm
A812610	A812612	Sm/Med	To 7.6
A812611	A812613	Lge/XL	7.6 +

For size, measure width of 2nd to 5th MCP.

## Rolyan® Rolyan progressive palm protector splint (with optional wrist support)

DVA RAP code AR03



Palm Splint

### Features:

- progressive soft splint for neurological conditions and arthritis
- adjust splint as condition improves by removing foam palm piece and replacing it with a larger series (see below)
- variable strapping system allows strapping thumb only, fingers only, or both thumb and fingers
- splint base is large enough to accommodate all hand sizes and can be modified for a smaller hand by trimming the edges
- wrist support style has a removable Ezeform stay that can be moulded for wrist extension and D-ring attachments for secure fit

- available in universal adult size
- each splint comes with one foam palm piece
- hand wash splint, air dry

### Application:

to determine which series to use-

- extend the fingers as far as comfortably tolerated
- measure flexion angle at MCP joint of tightest finger
- choose appropriate splint for current patient needs
- order by series/hand (eg series 1/left)



Palm Splint with Wrist Support

			Palm Splint		Palm Splint with Wrist Support		Replacement Pads
Series	Foam Size cm	MCP Flexion	Left	Right	Left	Right	
1	2.5	90	A812313	A812314	A812301	A812302	A812307
2	5.1	60	A812315	A812316	A812303	A812304	A812309
3	7.6	30	A812317	A812318	A812305	A812306	A812311

### Indications:

Series 1	Tight hands with moderate to severe spasticity or contractures
Series 2	Hands with mild to moderate spasticity or contractures
Series 3	Hands with mild spasticity or contractures, or non-spastic hands, to provide a functional position in a soft splint

## Rolyan® Rolyan Kwik-Form progressive hand splint

DVA RAP code AR03



### Features:

- can be custom fitted without tools or heat - simply bend the wire core to the desired position
- supports the hand and wrist and provides progressive extension of the fingers
- includes 3 hand rolls - as condition progresses, simply replace the hand roll
- for especially tight hands, can be used without hand roll
- removable cover sewn from comfortable, cool soft material for maximum compliance

- features hook and loop closures and removable machine washable cover
- includes laundry bag for washing and storing
- to determine size, measure wrist circumference
- bilateral design fits right or left

Splint	Replacement Covers	Size	Wrist
A812620	A812622R	Sm/Med	To 18
A812621	A812623R	Med/Lge	18cm +

Measure (cm) around wrist. Fits left or right

## Rolyan® Rolyan AirPro air graduate splint

DVA RAP code AR03



### Features:

- Air Therapy – low load prolonged stretch
- conforms to the shape of the hand accommodating different degrees of finger contracture
- ideal for severely contracted fingers

- removable cone shaped palmar air bladder can be used independently of wrist support
- extra cover included
- universal size

Left	Right
OC-AP31AGL	OC-AP31AGR

Universal size, left and right.

## Rolyan® Rolyan SoftPro Grip

DVA RAP code AR03



### Features:

- treats moderate to severe contractures of the wrist/hand/fingers
- easily moulded without tools yet sturdy enough to maintain shape during wear
- easily adjustable for ulnar and/or radial deviation

- finger separator and ulnar drift strap available separately to allow custom fitting
- extra cover included with every SoftPro splint
- palmar bar comfortably fits into the palm of a fist
- two incrementally larger palmar rolls included to gradually increase extension of the MCP and PIP joints

Left	Right	Size	Width	Length
OC-SP31FGSL	OC-SP31FGSR	Small	To 7.6	To 15.2
OC-SP31FGML	OC-SP31FGMR	Medium	7.6-8.9	15.2-20.3
OC-SP31FGLL	OC-SP31FGLR	Large	8.9-10.2	20.3-24.1
OC-SP31FGCL	OC-SP31FGCR	Child	5-7.6	17.8
OC-SP31FGYL	OC-SP31FGYR	Youth	5-7.6	20.3

Width: width of MP joints (cm). Length: wrist to finger tip (cm). Order by small/left etc.

## Rolyan® Rolyan Softpro finger separator (accessory to SoftPro splint)

DVA RAP code AR03



### Features:

- for use with SoftPro Grip and SoftPro functional resting hand splints
- provides finger separation and mild abduction
- protects finger skin integrity
- universal size

Product Code	Size
OC-SPFSS	Universal size

## Rolyan® Rolyan Softpro ulnar drift strap (accessory to SoftPro splint)

DVA RAP code AR03



### Features:

- for use with SoftPro Grip and SoftPro functional resting hand splints
- holds fingers in comfortable, neutral position
- enhances gentle three point leverage over the PIP joints
- universal size

Product Code	Size
OC-SPFSU	Universal size



## OPPO 4182 wrist/hand positioning splint

DVA RAP code AR03



### Features:

- two unstretchable straps give greater wrist stability
- incorporating a molded metal palmar stay which is used to position the wrist and hand
- insert each two fingers, except thumb, into the two unstretchable straps
- be sure the metal stay is on palm side of hand, then wrap the strap over hand and secure snugly
- thread through the buckles and secure the two hook and loop straps to desired compression

### Application:

- tendon repair
- post surgical use
- peripheral nerve injury

Left	Right	Size	Wrist (cm)
OPP4182S L	OPP4182S R	S	13.1-17.1
OPP4182M L	OPP4182M R	M	17.1-21.6
OPP4182L L	OPP4182L R	L	21.6-30.5

Measure wrist circumference.



## Rolyan Kwik-Form functional resting splint

DVA RAP code AR03



### Features:

- can be custom fitted without tools or heat - simply bend the wire core to the desired position
- supports the wrist, thumb and fingers in a functional or resting position

- removable cover sewn from comfortable, cool soft material for maximum compliance
- features hook and loop closures and removable machine washable cover
- includes laundry bag for washing and storing
- to determine size, measure wrist circumference

Splint		Replacement Covers			
Left	Right	Left	Right	Size	Wrist
A812701	A812700	A812705R	A812704R	Sm/Med	To 18
A812703	A812702	A812707R	A812706R	Med/Lge	18cm +

Measure (cm) around wrist.



## Rolyan SoftPro functional resting hand

DVA RAP code AR03



### Features:

- treats moderate contractures of the wrist/hand/fingers
- easily moulded without tools yet sturdy enough to maintain shape during wear
- easily adjustable for ulnar and/or radial deviation
- finger separator and ulnar drift strap available separately to allow custom fitting

- extra cover included with every SoftPro splint
- extra long wrist straps holds the wrist in the proper position and applies a gentle stretch
- thumb is "Bend to Fit" to support, position, and/or correct alignment of the thumb

Left	Right	Size	Width	Length
OC-SP31FFRSL	OC-SP31FFRSR	Small	To 7.6	To 15.2
OC-SP31FFRML	OC-SP31FFRMR	Medium	7.6-8.9	15.2-20.3
OC-SP31FFRLL	OC-SP31FFRLR	Large	8.9-10.2	20.3-24.1
OC-SP31FFRCL	OC-SP31FFRCR	Child	5-7.6	17.8
OC-SP31FFRYL	OC-SP31FFRYR	Youth	5-7.6	20.3

Width: width of MP joints (cm). Length: wrist to finger tip (cm)



## Cock-up resting hand splint

DVA RAP code AR03



### Features:

- pre-shaped plastic with velcro fastening straps
- use for carpal tunnel syndrome, RSI, resting splint
- available with or without thumb support
- s/m suitable for average female/small male hand
- large suitable for large female/average male hand

No Thumb			
Left	Right	Size	Width (cm)
B211LM	B211RM	S/M	7.5
B211LL	B211RL	Lge	9.0

With Thumb			
Left	Right	Size	Width (cm)
B212LM	B212RM	S/M	7.5
B212LL	B212RL	Lge	9.0