



Body Assist 250 deluxe wrist support – Basic support

Elastic

DVA RAP code

AR03



Features:

- firm and rigid anatomically shaped, compressive wrist brace
- adjustable thumb strap

Application:

- wrist sprain
- weak wrist

Product Code	Size	cm
NH250 1	S	To 16.5
NH250 2	M	16.5-19
NH250 3	L	19+

Measure wrist circumference.



McDavid A451 thermal wrist wrap – Basic support

Neoprene

DVA RAP code

AR03



Features:

- multi-purpose wrist wrap providing therapeutic warmth, firm compression and some protection
- thumb hole, fully adjustable
- velcro closure
- colour black
- one size fits all

Application:

- wrist sprain
- weak wrist

Product Code	size
NHA451	universal



OPPO 1081 wrist wrap – Basic support

Neoprene

DVA RAP code

AR03



Features:

- helps relieve pain through pressure on injured wrists
- adjustable hook and loop strap ensures comfortable fit and stabilization during activities
- provides firm support around weak or strained area

Application:

- wrist sprain
- weak wrist

Product Code	size	cm
OPP1081	universal	12-23

Measure wrist circumference.



OPPO 1083 wrist wrap w/thumb loop – Basic support

Neoprene

DVA RAP code

AR03



Features:

- provides comfortable compression and support for weak or stressed wrists
- helps treatment and prevention of injury
- retains warmth, which may aid and speed up recovery period

Application:

- wrist sprain
- weak wrist

Product Code	size	cm
OPP1083	universal	14-23

Measure wrist circumference.


OPPO 2081 wrist wrap – Basic support
DVA RAP code AR03

Features:

- helps relieve pain through pressure of injured wrists
- elastic hook and loop strap insures comfortable fit and stabilization during activities
- provides firm support around weak or strained area

Application:

- wrist sprain
- weak wrist

Product Code	size	cm
OPP2081	universal	13-23

Measure wrist circumference.


OPPO 2083 wrist wrap w/thumb loop – Basic support

Cotton blend

DVA RAP code AR03

Features:

- universal design fits left and right hand
- provides comfortable compression and support

Application:

- wrist sprain
- weak wrist

Product Code	size	cm
OPP2083	universal	14-23

Measure wrist circumference.


OPPO 2481 wrist wrap w/silicone – Basic support
DVA RAP code AR03

Features:

- provides localized support and compression, and avoids unnecessary compression to areas where it is not required
- tension is fully adjustable to provide from light support to total compression
- cohesive fabric, causing no pain or irritation to the skin whilst in use, comfortable to wear and painless to remove
- 7.5cm width
- secure the hook closure

Application:

- wrist sprain
- weak wrist

To Apply:

- place wrap on wrist and let the silicone side face to skin
- wrap around the wrist as close to the hand as possible
- apply firm pressure without impeding circulation

Product Code	size
OPP2481	universal


OPPO 1082 wrist splint – Stabilising

Neoprene

DVA RAP code AR03

Features:

- provides compression for weak or overstressed wrists
- spoon stay can be removed as healing progresses
- stay can be relocated to fit either left or right hand

Application:

- carpal tunnel syndrome
- post traumatic use

To Apply:

- position the contoured splint on the palm side of hand
- secure each hook and loop straps to desired compression
- the top strap goes between the thumb and forefinger and attaches at the back of hand

Product Code	Size	cm
OPP1082 1	S	14.0-15.9
OPP1082 2	M	15.9-18.4
OPP1082 3	L	18.4-21.0
OPP1082 4	XL	21.0-23.5

Measure wrist circumference.



OPPO 2288 wrist splint with strap

Cotton

DVA RAP code

AR03



Features:

- semi-rigid thumb spica portion is contoured to provide support, immobilization and thumb abduction
- provides wrist stability with moderate motion restriction
- an elastic wrist strap gives greater wrist stability
- aluminum stay can easily be reversed to fit right or left hand

Application:

- carpal tunnel syndrome
- post surgical use
- Colle's fracture

To Apply:

- position the contoured stay on the palm side of hand
- the top strap goes between the thumb and forefinger and attaches at the back of hand
- secure each hook and loop straps to desired compression
- wrap and secure the reinforced strap to desired compression

Product Code	Size	cm
OPP2288 1	S	14.0-15.9
OPP2288 2	M	15.9-18.4
OPP2288 3	L	18.4-21.0
OPP2288 4	XL	21.0-23.5

Measure wrist circumference.



OPPO 3082 wrist support with D-rings

Polyester

DVA RAP code

AR03



Features:

- hook and loop closures assure proper fit and easy application
- removable aluminum splint helps relieve pain.
- breathable lightweight fabric enhances comfort.
- elastic hook and loop strap may be used to limit wrist motion

Application:

- carpal tunnel syndrome
- Colle's fracture
- post surgical use

Left	Right	Size	cm
OPP3082 1 L	OPP3082 1 R	S	14.0-15.9
OPP3082 2 L	OPP3082 2 R	M	15.9-18.4
OPP3082 3 L	OPP3082 3 R	L	18.4-21.0
OPP3082 4 L	OPP3082 4 R	XL	21.0-23.5

Measure wrist circumference.



OPPO 2182 wrist splint with strap

Cotton

DVA RAP code

AR03



Features:

- an elastic wrist strap gives greater wrist stability
- reversible aluminum splint may be easily be relocated to fit right or left hand

Application:

- carpal tunnel syndrome
- severe wrist sprain

To Apply:

- position the contoured stay on the palm side of hand
- the top strap goes between the thumb and forefinger and attaches at the back of hand

- secure each hook and loop straps to desired compression
- wrap and secure the reinforced strap to desired compression

Product Code	Size	cm
OPP2182 1	S	14.0-15.9
OPP2182 2	M	15.9-18.4
OPP2182 3	L	18.4-21.0
OPP2182 4	XL	21.0-23.5

Measure wrist circumference.



Body Assist 260 deluxe wrist/hand brace

DVA RAP code

AR03



Features:

- contoured support for weak, injured wrists
- pliable rigid splint to minimize wrist movement

Application:

- for fractures
- after removal of plaster cast
- injuries to ligaments, muscles and tendons

Product Code	Size	cm
NH260 1	S	To 16.5
NH260 2	M	16.5-19
NH260 3	L	19-21.5
NH260 4	XL	22+

Measure wrist circumference.



McDavid A454 thermal carpal tunnel wrist splint

Neoprene

DVA RAP code **AR03**



Features:

- contoured thermal neoprene for weak injured wrists
- includes a plantar splint to minimize wrist movement
- specify left or right
- one size fits all

Application:

- for use with fractures, injuries to muscles, ligaments and tendons

Product Code		Size
NHA454L	Left	universal
NHA454R	Right	universal



Rehband 4024 manu rheuma flex wrist

DVA RAP code **AR03**



Features:

- wrist support made of soft, skincaring and breathable material
- supports the wrist through a flexible palmar splint
- velcro closure simplifies applying the support for rheumatic patients
- a soft and pre-shaped plastic splint provides light support
- easy application through D-rings
- straight model to fit the arthritic hand
- the placement of the velcro is chosen to avoid any pressure on caput ulnae

Application:

- rheumatoid arthritis
- osteoarthritis

Left	Right	Size	cm
REH4024L XS	REH4024R XS	XS	14-15
REH4024L S	REH4024R S	S	15-16.5
REH4024L M	REH4024R M	M	16.5-18
REH4024L L	REH4024R L	L	18-19.5
REH4024L XL	REH4024R XL	XL	19.5-21
REH4024L XXL	REH4024R XXL	XXL	21-22.5

Measure wrist circumference.



Rehband 4055 manu comfort T Wrist

DVA RAP code **AR03**



Features:

- stabilising wrist brace
- easily applied
- excellent comfort through the temperature regulating features of ComforTemp
- includes plastic and steel splints
- 3cm shorter than REH4058

Application:

- artralgia
- tenosynovitis
- post-fracture management

Left	Right	Size	cm
REH4055L 1	REH4055R 1	XS-S	13-16
REH4055L 2	REH4055R 2	M-L	16-19
REH4055L 3	REH4055R 3	XL-XXL	19-23

Measure wrist circumference.



Rehband 4058 manu comfort T stable wrist

DVA RAP code **AR03**



Features:

- extra stabilising wrist brace
- easily applied
- excellent comfort through the temperature regulating features of ComforTemp
- includes plastic and steel splints
- 3cm longer than REH4055

Application:

- post-trauma/post-operative management
- partial paralysis

Left	Right	Size	cm
REH4058L 1	REH4058R 1	XS-S	13-16
REH4058L 2	REH4058R 2	M-L	16-19
REH4058L 3	REH4058R 3	XL-XXL	19-23

Measure wrist circumference.



Rehband 4062 manuell EQ wrist support

DVA RAP code AR03



Features:

- long-lasting and dirt-repelling material
- plastic splint
- closed thumb grip
- pain relieving and unloading effects
- may be used with #401 stabilising strap

Application:

- used to stabilise and support the over-used wrist

Left	Right	Size	cm
REH4062L XS	REH4062R XS	XS	14-15
REH4062L S	REH4062R S	S	15-16.5
REH4062L M	REH4062R M	M	16.5-18
REH4062L L	REH4062R L	L	18-19.5
REH4062L XL	REH4062R XL	XL	19.5-21
REH4062L XXL	REH4062R XXL	XXL	21-22.5

Measure wrist circumference.



Rehband 4063 rheuma manuell EQ wrist support

DVA RAP code AR03



Features:

- stabilising wrist brace
- durable and dirt-repelling material
- adjustable thumb grip
- includes plastic splint

Application:

- rheumatoid arthritis
- arthralgia
- over-use syndrome

Left	Right	Size	cm
REH4063L XS	REH4063R XS	XS	14-15
REH4063L S	REH4063R S	S	15-16.5
REH4063L M	REH4063R M	M	16.5-18
REH4063L L	REH4063R L	L	18-19.5
REH4063L XL	REH4063R XL	XL	19.5-21
REH4063L XXL	REH4063R XXL	XXL	21-22.5

Measure wrist circumference.



Rehband 4074 diagonal lycra short wrist

DVA RAP code AR03



Features:

- a thumb support! with splint to stabilize the CMC- and MCP-joints of the thumb
- the material offers optimal fit and comfort
- pain relieving and unloading effects

Application:

- used for rheumatic or degenerative problems

Left	Right	Size	cm
REH4074L XS	REH4074R XS	XS	13-14
REH4074L S	REH4074R S	S	14-16
REH4074L M	REH4074R M	M	16-18
REH4074L L	REH4074R L	L	18-20
REH4074L XL	REH4074R XL	XL	20-22

Measure wrist circumference.



Rehband 4088 wristform – Heat mouldable

DVA RAP code AR03



Features:

- provides exceptional stability while immobilizing the wrist
- made of thermoprene material for optimal comfort
- a firm wrist stay supports the palmar side, while a dorsal pocket holds a thermoplastic insert
- both the stay and insert can be heat-moulded

Application:

- ideal for post-scapoid fractures, tendinitis, epicondylitis, carpal tunnel syndrome and general instabilities or hyperextension problems of the wrist

To Apply:

- pour hot water into the pocket to soften material
- pour water out, blot dry and mould

Left	Right	Size	cm
REH4088L S	REH4088R S	S	14-16
REH4088L M	REH4088R M	M	16-18
REH4088L L	REH4088R L	L	18-20
REH4088L XL	REH4088R XL	XL	20-22

Measure wrist circumference.


Rehband 4142 manu 3D

Elastic

DVA RAP code AR03

Features:

- high quality 3 tab wrist brace with flexion stay
- provides a light support to the basal thumb joint

Application:

- various thumb and wrist disorders

Left	Right	Size	cm
REH4142L XS	REH4142R XS	XS	14.5-15.5
REH4142L S	REH4142R S	S	15.5-16
REH4142L M	REH4142R M	M	16-17
REH4142L L	REH4142R L	L	17-18
REH4142L XL	REH4142R XL	XL	18-19
REH4142L XXL	REH4142R XXL	XXL	19-20

Measure wrist circumference.


Rehband 4143 manu 3D basis

Elastic

DVA RAP code AR03

Features:

- high quality 3 tab wrist brace with flexion stay
- allows free thumb motion
- six sizes left and right

Application:

- various thumb and wrist disorders

Left	Right	Size	cm
REH4143L XS	REH4143R XS	XS	14.5-15.5
REH4143L S	REH4143R S	S	15.5-16
REH4143L M	REH4143R M	M	16-17
REH4143L L	REH4143R L	L	17-18
REH4143L XL	REH4143R XL	XL	18-19
REH4143L XXL	REH4143R XXL	XXL	19-20

Measure wrist circumference.


Rolyan D-ring wrist brace
DVA RAP code AR03

Features:

- non-elastic construction does not restrict circulation
- circumferential design holds the brace in position during application
- D-ring closures provide tension control to adjust amount of support
- double layer of cotton stockinette under closures cushions straps
- forearm length is proportional to brace size
- removable, adjustable aluminium support bar
- splints available in regular (A611) and long (A915) length
- regular length splint is also available with a larger thumb hole, which accommodates larger thumbs and fluctuating oedema.
- fabricated from a polyester/cotton fabric laminated to 1cm foam

Application:

- for high-compliance treatment of arthritis, carpal tunnel syndrome, tendonitis and other wrist disorders

Left	Right	Size	cm	Length
Regular Fit				
A611LXS	A611RXS	XS	To 14.6	15.9
A611LS	A611RS	S	14.6-16.5	16.5
A611LM	A611RM	M	16.6-19.1	17.8
A611LL	A611RL	L	19.2-21.6	19.0
A611LXL	A611RXL	XL	Over 21.6	19.7
Regular Fit - Larger thumb opening				
A611101	A611100	XS	To 14.6	15.9
A611103	A611102	S	14.6-16.5	16.5
A611105	A611104	M	16.6-19.1	17.8
A611107	A611106	L	19.2-21.6	19.0
A611109	A611108	XL	Over 21.6	19.7
Long Fit				
A915LXS	A915RXS	XS	To 14.6	21.0
A915LS	A915RS	S	14.6-16.5	21.6
A915LM	A915RM	M	16.6-19.1	22.2
A915LL	A915RL	L	19.2-21.6	22.6
A915LXL	A915RXL	XL	Over 21.6	23.2

Measure wrist circumference.



Rolyan D-ring wrist long fit - firm

DVA RAP code **AR03**



Features:

- volar and dorsal stays for extra wrist support
- prevents unwanted flexion or extension, yet allows unobstructed motion of the fingers and thumb
- non-elastic circumferential design for easy application and for firm, evenly distributed pressure without restricting circulation
- protective padding lined with cotton stockinette
- D-ring closure system for simple application, removal and customized tension control
- double layer of stockinette under closures cushions straps
- forearm length is proportional to brace size
- removable, adjustable aluminum bar for additional support
- machine wash and air dry

Application:

- for high-compliance treatment of arthritis, carpal tunnel syndrome, tendonitis and other wrist disorders

Left	Right	Size	cm
A600301	A600300	XS	To 14.6
A600303	A600302	S	14.6-16.5
A600305	A600304	M	16.6-19.1
A600307	A600306	L	19.2-21.6

Measure wrist circumference.



Rolyan wrist support

DVA RAP code **AR03**



Features:

- creep resistant elastic which prevents bunching that can restrict circulation
- anatomical design provides superior comfort and function and the two piece design assures optimal fit
- includes a removable adjustable stay
- hand or machine wash, air dry

Application:

- for painful wrists or to prevent unwanted wrist flexion or extension

Left	Right	Size	cm
A619063	A619062	S	14.0-16.5
A619065	A619064	M	16.6-18.4
A619067	A619066	L	18.5-20.3
A619069	A619068	XL	20.4-22.2

Measure wrist circumference.



Rolyan Takeoff universal wrist splint

DVA RAP code **AR03**



Features:

- semi-rigid foam support with plastic stay
- soft, breathable Fabrifoam material provides friction to keep splint from migrating and helps remove moisture
- one size fits most adults

Application:

- mild to moderate wrist sprains and strains, carpal tunnel syndrome

Left	Right	Size
A919801	A919800	universal



Rolyan Alignrite wrist support w/wrist strap

DVA RAP code **AR03**



Features:

- provides excellent support for the wrist while allowing fingers and thumb free movement
- woven cotton-coated rubber fabric has open weave for ventilation
- wrap-around wrist strap provides extra support
- hand or machine wash, air dry
- latex free

Application:

- helps reduce pain and inflammation from conditions such as arthritis, sprains, and cumulative trauma injuries such as carpal tunnel syndrome

Left	Right	Size	cm
A619001	A619002	XS	11.4-14.0
A619003	A619004	S	14.0-16.5
A619005	A619006	M	16.6-18.4
A619007	A619008	L	18.5-20.3
A619009	A619010	XL	20.4-22.2

Measure wrist circumference.


Rolyan wrist splint – Firm
DVA RAP code AR03

Features:

- unlike neoprene or perforated neoprene, breathoprene allows air through, and wicks perspiration away from the skin for enhanced comfort and compliance
- the two-way stretch 3.2mm thick, soft material provides compression and a perfect fit
- removable metal dorsal and palmar stays and wrist D-Ring strap for added support
- wraparound design
- latex free

Application:

- helps reduce pain and inflammation from conditions such as arthritis, sprains, and cumulative trauma injuries such as carpal tunnel syndrome

Left	Right	Size	cm
DM45101	DM45102	S	14.0-16.5
DM45103	DM45104	M	16.6-18.4
DM45105	DM45106	L	18.5-20.3
DM45107	DM45108	XL	20.4-22.2

Measure wrist circumference.


OPPO 1289 wrist/thumb support w/strap

Coolprene

DVA RAP code AR03

Features:

- three removable metal splints in thumb area for improved stability
- additional strap for stabilization during activities
- supports the wrist and thumb, leaving the palm free
- ideal for sprains of thumb, ligament support, and minor fractures of the thumb

Application:

- De Quervain's disease
- thumb sprain/strain
- CMC joint injury

To Apply:

- position thumb on the thumb locking area of support
- be sure the metal splints are on the lateral side of thumb
- secure the hook and loop strap beside thumb area to desired compression
- pull the strap around the wrist and secure hook and loop strap to desired compression

Product Code	Size	cm
OPP1289 1	S	13.0-15.2
OPP1289 2	M	15.2-17.2
OPP1289 3	L	17.2-19.1
OPP1289 4	XL	19.1-21.0

Measure wrist circumference.


OPPO 2184 wrist/ thumb sleeve

Cotton blend

DVA RAP code AR03

Features:

- pull-on sleeve in four-way stretch tubular elastic
- not splinted, light compression
- universal design fits left and right hand

Application:

- weak thumb and wrist
- wrist sprain/ strain

Product Code	Size	cm
OPP2184 1	S	14.0-15.9
OPP2184 2	M	15.9-18.4
OPP2184 3	L	18.4-21.0
OPP2184 4	XL	21.0-23.5

Measure wrist circumference.


OPPO 3182 wrist support with thumb

Polyester

DVA RAP code AR03

Features:

- wraparound design and hook and loop straps for easy fit.
- supports and limits the motion of the thumb and wrist joint
- elastic fabric provides compression.
- elastic hook and loop strap may be used to limit wrist motion

Application:

- carpal tunnel syndrome
- Colle's fracture
- post surgical use

Left	Right	Size	cm
OPP3182 1 L	OPP3182 1 R	S	14.0-15.9
OPP3182 2 L	OPP3182 2 R	M	15.9-18.4
OPP3182 3 L	OPP3182 3 R	L	18.4-21.0
OPP3182 4 L	OPP3182 4 R	XL	21.0-23.5

Measure wrist circumference.



OPPO 1089 wrist thumb support

Neoprene



Features:

- thumb orthosis especially suitable for rheumatics
- relieves pain and stabilizes CMC joint
- reinforced metal stick stabilizes thumb
- adjustable hook and loop fasteners insure comfortable fit and stabilization during activities
- thermal wrist support for tenosynovitis, carpal tunnel syndrome and general wrist support

Application:

- De Quervain's disease
- thumb sprain/strain
- CMC joint injury

To Apply:

- position thumb on the thumb locking area of support
- be sure the metal splint is on the lateral side of thumb
- secure the hook and loop strap beside thumb area to desired compression
- secure elastic hook and loop fastener to desired compression

Product Code	Size	cm
OPP1089 1	S	14.0-15.9
OPP1089 2	M	15.9-18.4
OPP1089 3	L	18.4-21.0
OPP1089 4	XL	21.0-23.5

Measure wrist circumference.

DVA RAP code AR03



OPPO 1188 wrist/thumb support – 3 lengths - Heat mouldable

Neoprene



Features:

- relieves pain and stabilizes thumb
- thermoplastic support stabilizes thumb
- thermoplastic insert may be removed as healing progresses
- keeps the CMC and MCP joints in a neutral position
- available in 3 lengths (14cm, 16.5cm (pictured) and 21cm)

Application:

- ulnar/radial collateral ligament injury
- arthritis
- DeQuervain's tenosynovitis

- MCP or CMC joint dislocation
- gamekeeper's thumb

To Apply:

- insert the thermoplastic plate into the pocket
- pour hot water (about 75) into the pocket
Water temperature cannot be above 100
- pour out water when thermoplastic plate is mouldable
- mould the plate and wait until it becomes stiff
- dry support with a towel
- adjust the Velcro strap for correct fitting of the support

Short 14cm		Medium 16.5cm		Long 21cm		Size	cm
Left	Right	Left	Right	Left	Right		
OPP1188-5 1 L	OPP1188-5 1 R	OPP1188-6 1 L	OPP1188-6 1 R	OPP1188-8 1 L	OPP1188-8 1 R	S	14.0-15.9
OPP1188-5 2 L	OPP1188-5 2 R	OPP1188-6 2 L	OPP1188-6 2 R	OPP1188-8 2 L	OPP1188-8 2 R	M	15.9-18.4
OPP1188-5 3 L	OPP1188-5 3 R	OPP1188-6 3 L	OPP1188-6 3 R	OPP1188-8 3 L	OPP1188-8 3 R	L	18.4-21.0
OPP1188-5 4 L	OPP1188-5 4 R	OPP1188-6 4 L	OPP1188-6 4 R	OPP1188-8 4 L	OPP1188-8 4 R	XL	21.0-23.5
OPP1188-5 5 L	OPP1188-5 5 R	OPP1188-6 5 L	OPP1188-6 5 R	OPP1188-8 5 L	OPP1188-8 5 R	XXL	23.5-26.0

Measure wrist circumference.

DVA RAP code AR03



Rehband 4025 rheuma wrist & thumb support



Features:

- wrist/thumb support made of soft, skincaring and breathable material
- supports the wrist through a flexible palmar splint
- supports the thumb by aluminium stay
- velcro closure simplifies applying the support for rheumatic patients
- easy application through D-rings
- straight model to fit the arthritic hand
- the placement of the velcro is chosen to avoid any pressure on caput ulnae

Application:

- rheumatoid arthritis
- osteoarthritis

Left	Right	Size	Cm
REH4025L XS	REH4025R XS	XS	14-15
REH4025L S	REH4025R S	S	15-16.5
REH4025L M	REH4025R M	M	16.5-18
REH4025L L	REH4025R L	L	18-19.5
REH4025L XL	REH4025R XL	XL	19.5-21
REH4025L XXL	REH4025R XXL	XXL	21-22.5

Measure wrist circumference.

DVA RAP code AR03



Rehband 4146 manu 3D pollex



Features:

- high quality 3 tab wrist brace with flexion stay
- aluminium thumb stay
- six sizes left and right

Application:

- various thumb and wrist disorders

Left	Right	Size	cm
REH4146L XS	REH4146R XS	XS	14.5-15.5
REH4146L S	REH4146R S	S	15.5-16
REH4146L M	REH4146R M	M	16-17
REH4146L L	REH4146R L	L	17-18
REH4146L XL	REH4146R XL	XL	18-19
REH4146L XXL	REH4146R XXL	XXL	19-20

Measure wrist circumference.

DVA RAP code AR03



Rolyan D-ring wrist/thumb spica

DVA RAP code AR03

Features:

- made of cotton/polyester fabric lined with cotton stockinette
- circumferential design with D-ring closures makes the brace easy to put on and comfortable to wear
- two metal stays support wrist and thumb
- adjust the amount of wrist extension with aluminium stay
- latex free

Application:

- allows healing after injury and decreases pain and inflammation of painful wrists and thumb tendonitis

Left	Right	Size	cm
A6151	A6152	XS	To 14.6
A6153	A6154	S	14.6-16.5
A6155	A6156	M	16.6-19.1
A6157	A6158	L	19.2-21.6
A6159	A61510	XL	Over 21.6

Measure wrist circumference.



Rolyan D-ring wrist/thumb spica - thermal

DVA RAP code AR03

Features:

- made from a neoprene blend with a Trioxon lining
- maintains body heat, yet allows skin to breathe
- circumferential design with D-ring closures makes the brace easy to put on and comfortable to wear
- two metal stays support wrist and thumb
- adjust the amount of wrist extension with aluminium stay
- latex free

Application:

- allows healing after injury and decreases pain and inflammation of painful wrists and thumb tendonitis

Left	Right	Size	cm
A61511	A61512	XS	To 14.6
A61513	A61514	S	14.6-16.5
A61515	A61516	M	16.6-19.1
A61517	A61518	L	19.2-21.6
A61519	A615120	XL	Over 21.6

Measure wrist circumference.



Body Assist N29 thumb abduction splint

Coolprene

DVA RAP code AR03

Features:

- minimises thumb movement while leaving palm and fingers free

Application:

- used in treatment of strains and sprains to metacarpal joint

Product Code	Size	cm
NHN29 1	S/M	To 19
NHN29 3	L/XL	19 +

Measure wrist circumference.



Body Assist 280 thumb splint

DVA RAP code AR03

Features:

- designed to minimize wrist and thumb movement
- supplied with two interchangeable splints – a rigid splint, and pliable metal stay

Application:

- useful in the treatment of strains and sprains of the thumb's metacarpal joint

Product Code	Size	cm
NH280 1	S	To 16.5
NH280 2	M	16.5-19
NH280 3	L	19-21.5
NH280 4	XL	22+

Measure wrist circumference.



Rehband 4064 EQ diagonal lycra thumb

DVA RAP code AR03



Features:

- a thumb support made of durable lycra
- has a splint to stabilise the CMC - and MCP - joints of the thumb
- the elasticity of the material is used diagonally for optimal fit and support

Application:

- provides pain relief and unloading to the thumb

Left	Right	Size	cm
REH4064L XS	REH4064R XS	XS	13-14
REH4064L S	REH4064R S	S	14-16
REH4064L M	REH4064R M	M	16-18
REH4064L L	REH4064R L	L	18-20
REH4064L XL	REH4064R XL	XL	20-22

Measure wrist circumference.



Rehband thumboform splint – 3 lengths – Heat mouldable

DVA RAP code AR03



Note: long 4085

Features:

- the thumboform thumb support comes with a heat-mouldable insert
- this splint stabilises the CMC and MCP joint of the thumb
- materials consist of neoprene, polyamide jersey
- available in short (4084), long (4085) and x/long (4086) styles

Application:

- indicated for rheumatoid arthritis, osteoarthritis, or for support after injury

To Apply:

- place the plastic insert in the pocket
- pour hot water (approximately 75°C) into the pocket, wait until the plastic is mouldable, dry with a towel.
- apply the brace and close the velcro carefully.
- mould the plastic and wait until it gets stiff
- let dry and the brace is ready to use.
- apply the brace, close the velcro and adjust the thumb grip

Short		Long		X- Long			
Left	Right	Left	Right	Left	Right	Size	Cm
REH4084L XS	REH4084R XS	REH4085L XS	REH4085R XS	REH4086L XS	REH4086R XS	XS	13-14
REH4084L S	REH4084R S	REH4085L S	REH4085R S	REH4086L S	REH4086R S	S	14-16
REH4084L M	REH4084R M	REH4085L M	REH4085R M	REH4086L M	REH4086R M	M	16-18
REH4084L L	REH4084R L	REH4085L L	REH4085R L	REH4086L L	REH4086R L	L	18-20
REH4084L XL	REH4084R XL	REH4085L XL	REH4085R XL	REH4086L XL	REH4086R XL	XL	20-22

Measure wrist circumference.



Rolyan Collum CMC thumb brace

DVA RAP code AR03



Features:

- made of quality leather with a soft suede lining for comfort and durability
- restricts extremes of motion in the thumb joint, while allowing the user to perform routine hand functions
- the support fits over the thumb and wraps around the wrist
- it is held in place with two hook and loop fasteners
- the brace is available for the right or left hand

Application:

- designed for people with arthritis at the CMC and MCP joints of the thumb

Left	Right	Size	cm
A920LXS	A920RXS	XS	To 9.5
A920LS	A920RS	S	9.8-10.8
A920LM	A920RM	M	11.1-12.1
A920LL	A920RL	L	over12.4

Measure length from the top thumb-joint, along outside of thumb to the wrist-bone.

Rolyan® Rolyan gel thumb spica splint

DVA RAP code AR03


Features:

- cushions and helps protect the thumb.
- gel shell pad provides gentle pressure on the incision site to help reduce scar formation and hypersensitivity
- optional polyform stay along radial side of thumb provides additional support to the thumb
- gel pad cushions area over first dorsal compartment to help protect it and disperse energy of impacts
- hand wash, air dry

Application:

- provides a barrier against trauma, and semi-rigid support to CMC and MCP joints

Left	Right	Size	cm
A599001	A599002	XS	To 9.5
A599003	A599004	S	9.5-10.8
A599005	A599006	M	10.9-12
A599007	A599008	L	over12.1

Measure length from the top thumb-joint, along outside of thumb to the wrist-bone.

Rolyan® Rolyan Takeoff thumb splint

DVA RAP code AR03


Features:

- made using fabriform NuStimWrap
- the splint is available with either a standard plastic thumb stay for increased support, or a mouldable stay for custom fit
- provides warmth, compression and support to reduce pain
- easy to apply over joints and helps prevent overuse
- One size fits most

Application:

- recommended for partial immobilisation of the thumb after surgery or trauma
- the mouldable stay is recommended for DeQuervain's tenosynovitis, skier's thumb, gamekeeper's thumb, collateral ligament support, avulsion fractures of the thumb, and arthritis

Style	Left	Right
Regular	A952713	A952712
With Mouldable Insert	A952701	A952700

Rolyan® Rolyan thermo-form thumb splint – Heat mouldable

DVA RAP code AR03


Features:

- these splints contain a pocket on the inside where the thermoplastic piece is inserted.
- the thermoplastic material can be moulded directly in the pocket by pouring in hot water
- the thermoplastic material becomes mouldable in approx. 1 minute, allowing the therapist to mould the splint to the thumb for proper positioning
- the splint comes with one strap

Application:

- thumb arthritis (osteo or rheumatoid)
- MCP joint injuries (ligament, arthritis, and tendon lacerations)
- MCP fusion, tendonitis, joint fractures

Left	Right	Size	cm
C55006101	C55006102	XS	To 14
C55006103	C55006104	S	14-16.5
C55006105	C55006106	M	16.5-19.1
C55006107	C55006108	L	19.1-21.6
C55006109	C550061010	XL	21.6-23

Measure wrist circumference.

Rolyan® Rolyan Action thumsling

DVA RAP code AR03


Features:

- four-way stretch material conforms in the palm without adding bulk
- reduces CMC subluxation and pain without limiting wrist or thumb motion
- the sling strap supports the thumb and assists in maintaining web space to reduce CMC subluxation

Application:

- support for CMC joints

Left	Right	Size	cm
3P-P2006-L23	3P-P2006-R23	S/M	15.0-18.4
3P-P2006-L34	3P-P2006-R34	M/L	18.5-23.5

Measure circumference at MP's.