

Table of Contents

- 2 Contact Information
- 3 About Swell Spots
- 4 Curved & Dorsal
- 5 Palmar & Oval
- 6 Katina & Sobol
- 7 Joju & Kimbe
- 8 Breast & Versi-Face
- 9 Thora-Thigh & Reynolds
- 10 Crescent & We'a Toe
- 11 Labia & Scrotum
- 12 Lateral Bra & Full Bra
- 13 LT Knee Shaper & Ring
- 14 Sweet Spot & Tic-Tac
- 15 Care Instructions & Ordering Information

Customer Service

Customer Service Representatives are available to assist you from 7 AM to 7 PM Central, Monday through Friday.

Tel: (855) 892-4140 *Fax*: (414) 892-4150

Email: info@solarismed.com

Solaris 3880 W Wheelhouse Road Milwaukee, WI 53208 USA

Swell Spots

Swell Spots are professionally designed treatment pads that tuck into garments or under bandages. Use them any time, day or night, to address particularly troublesome areas of edema, fill concavities, or diffuse pressure points. They can also be customized to provide a unique solution*.

Highlights



Proprietary fabric is comfortable, breathable, wickable, and thermal release.




Grey microfiber fabric grips skin and encourages gentle tissue stretch to help open the initial lymphatics and encourage reabsorbtion.



Help break up tough fibrotic tissue while protecting delicate skin.




Evenly disperse irregular pressures to prevent compression hot-spots.



- Dorsal SP-02
- 1) LT Knee Shaper SP-17
- ② Ring with Velcro SP-18-A
- Thora-Thigh SP-11
- 4 Oval SP-04
- Crescent SP-13


^{*} Swell Spots can be custom made. If your custom Swell Spot is put into standard production, you'll be given the option to name it before it goes into the catalog.



Curved SP-01

Primarily used for hard to accommodate areas (like the ankles and patella), but also effective around the elbow, for inferior breast scarring, around the eyes, and behind the ear.

Item Code	Width	Height
SP-01	5cm	13cm



Dorsal SP-02


Designed for use on the dorsum of either the hand or foot.

Item Code	Width	Height
SP-02-S	10cm	14cm
SP-02-L	10cm	19cm

Palmar SP-03

Primarily for use on palmar surface of hand; also effective on the dorsum, concave shape of palm, and lumpectomy scarring when used in a bra.


Item Code	Width	Height
SP-03-S	7cm	8cm
SP-03-L	8cm	9cm




Oval SP-04

Versatile and easy to use with a simple, multi-purpose design that can be placed in many areas, such as inside the hand splints, along the ulnar aspect of the forearm, and the groin.

Item Code	Width	Height
SP-04	8cm	19cm





Katina SP-05



Tapered shape useful in treating difficult areas, such as the ulnar aspect of the forearm, the IT band of the lateral thigh, and the anterior surface of the calf.

Item Code	Width	Height
SP-05-S	11cm	36cm
SP-05-L	18cm	42cm



x 18 cm



Sobol SP-06

Maria Sobo

Addresses post-mastectomy fibrotic tissue on chest wall due to radiation and/or surgical incision scarring. Also great for use on the hip, anterior or posterior knee, and abdomen when used under bandages and compression garments.

Item Code	Diameter
SP-06-S	18cm
SP-06-L	22cm

JoJu SP-07

Specifically addresses swelling of the thenar web space; Therapists have also used it to soften scars from carpal tunnel surgery.




Item Code	Width	Height
SP-07	4.5cm	7cm

esigned by

Kimbe SP-08

Appropriate for reducing swelling of the scapula, as well as over the groin and neck area.

Item Code	Width	Height
SP-08	14cm	24cm








Designed specifically for swelling and fibrosis of the breast area after a lumpectomy or reconstructive surgery. New seamless design allows SP-09 to be used on both sides of the chest without irritation.

Item Code	Cup Size*
SP-09-S	A – B cup
SP-09-M	C – D cup
SP-09-L	D – DD cup

^{*} Sizes are approximate



Versi-Face SP-10


Designed to address swelling around the eyes; also effective when placed under the chin covering the neck or rotated to cover the mandible.

Item Code	Width	Height
SP-10	29cm	16cm

Thora-Thigh SP-11

Designed to provide coverage for the chest, groin, medial thigh, and ankle areas; Helps address radiation fibrosis and disperse distal pressure of bandages and leg garments.

Item Code	Width	Height
SP-11	52cm	22cm



Reynolds SP-12


Reynolds strap is constructed with Velcro-receptive material that hooks to the Swell Spot; it can be placed in multiple locations around the head and neck.





Item Code	Spot Width	Spot Height
SP-12	16.5cm	7.5cm


51 cm



Crescent SP-13

Commonly used in skin folds; helps soften upper and lower lobules and supports lax abdominal tissue.

Item Code	Length	Height	Diameter
SP-13-S	51cm	6.5cm	9cm
SP-13-M	61cm	6.5cm	11cm
SP-13-L	71cm	6.5cm	12cm



We'a Toe SP-14


Exclusive digit spacers offer comfort and cradle individual toes. Clients must be able to reach toes to don.

Item Code	MTP Circumference
SP-14-S	22cm
SP-14-M	25cm
SP-14-L	32cm

Labia SP-15-F

Tapered design enhances comfort while applying lower abdominal coverage to help direct fluid away from genital area; can be worn inside or over underwear.

Item Code	Width	Height
SP-15-F	33cm	53cm







33 cm **X**

Scrotum SP-15-M

Designed to gently cradle scrotum; accommodates up to 27 cm of scrotal circumference (further customization available).

Item Code	Width	Height
SP-15-M	34cm	53cm









Designed to tuck inside compression bra pockets or any regular bra to reduce scarring and protect delicate areas from bra chest strap pressures.

Item Code	Width	Height
SP-16-A	10cm	16cm





Full Bra SP-16-B


Designed for use under a Belisse® bra or inside sports bras with front opening.

Item Code	Cup Size
SP-16-B-S	А – В сир
SP-16-B-M	B – C cup
SP-16-B-L	C – D cup

LT Knee Shaper SP-17

Use to pad bony prominences, diffuse bandaging pressures, and treat tissue induration on the knees or elbows. Helps contour/shape knee while protecting delicate structures of the popliteal fossa.

Item Code	Girth Range
SP-17	35 – 65cm



Ring SP-18

Designed to fill crevices or lift lobules at the ankle and wrist; may also be used for abdominal pannuses. Available with or without a Velcro® closure.

Item Code	Length	Height	Diameter
SP-18	60cm	6.5cm	9cm




60 cm





Addresses residual and transient pelvic floor swelling in pregnant and postpartum clients, as well as perineum and labia swelling. One-sided non-skid surface prevents shifting and can be worn inside, over, or between underwear.

Item Code	Width	Height
SP-19	7cm	15cm



SP-20-M

Tic-Tac SP-20

Works to address larger areas of fibrosis, such as the back, abdomen, or thigh.

Item Code	Size
SP-20-S	20 x 20cm
SP-20-M	30 x 30cm
SP-20-L	40 x 40cm

30

x 30 cm

How to Order

Solaris places a high importance on providing products exclusively through medical professionals. For a complete list of retail locations in your area, visit **solarismed.com**

If your facility provides durable medical equipment, call or e-mail to set up an account. Solaris Customer Service Representatives are available to assist you from 7 AM until 7 PM Central, Monday through Friday.

Return Policy & Warranty

Swell Spots may be returned if found to be defective within two weeks of delivery. Shipping fees are nonrefundable.

Swell Spots are not warranted.

Care Instructions

Washing

- Machine wash with mild detergent on permanent press.
- Double rinse to remove soap residue.
- Do not use fabric softener, bleach, or other additives.

Drying

Tumble dry on permanent press

Indications

Edema

Contraindications

- Arterial insufficiency or degeneration
- Deep Vein Thrombosis
- Untreated congestive heart failure
- · Localized or systemic infection
- Untreated cancer
- Absent or severely impaired sensation

Solaris

3880 W Wheelhouse Road Milwaukee, WI 53208 USA

Tel: (855) 892-4140 *Fax:* (414) 892-4150

Email: info@solarismed.com

Web: solarismed.com

- f /solarisinc
- @solarisinc
- /solarisinc